

LABEE CHECKIT GO!



The LabBee Checkit Go is a simple, stand-alone, disposable cartridge that allows you to check the accuracy of your liquid handling robot or multi-channel pipette in just seconds. No additional equipment is required.

BENEFITS:

SIMPLE, USER-FRIENDLY OPERATION THAT REQUIRES NO TRAINING

- Dispense the dye into the wells on the flip side of the cartridge.
- Flip the tab (or have the robot do it)
- Check and compare the level of the dye against the graduations

FAST AND ACCURATE

- Measure dispensed volumes of up to 8 channels simultaneously in under 10 seconds
- Directly measures the volume dispensed compared to other indirect methods (such as the colorimetric and gravimetric methods)
- Measurements are guaranteed to be accurate within 2%

WORKS WITH YOUR MULTI-CHANNEL PIPETTE AND YOUR AUTOMATED LIQUID HANDLER

- Compatible with all brands of multi-channel pipettes and electronic pipettes
- Compatible with all liquid handlers that handle well plates

NO CAPITAL INVESTMENT NEEDED

- No annual maintenance costs
- Fraction of the cost of unnecessary calibration
- No staff training hours required

VERSATILITY

- Facilitates the measurement of solutions with differing viscosities
 - Up to 10% DMSO
 - Up to 10% Ethanol
 - Up to 10% Glycerol
 - 100% Serum
 - 100% Plasma

ORDERING INFORMATION

SKU	Description
MSLO-8V5-6	5ul Checkit Go - box of 6
MSLO-8J10-6	10ul Checkit Go - box of 6
MSLO-8K20-6	20ul Checkit Go - box of 6
MSLO-8C50-6	50ul Checkit Go - box of 6

