



# Universal Manual & Robotic Tips

- *ISO-Certified Manufacturing*
- *Precise Results*
- *Always Available*



LABBEE is a new brand name built on many years' experiences in the life science industry meeting customer supply chain needs. Their mission is based around five core principles that outline their commitment guarantee to customers:

### **RELIABLE. AVAILABLE. TRUSTWORTHY. ECONOMICAL. DELIVERABLE.**

We aim to provide our partners with the product they need without any compromise on service or quality. For our first product range, we identified pipette tips as being one of the core needs in every laboratory. Our range covers both universal manual tips and robotic tips where market choice and price were limiting for customers until now. Our range of pipette tips are manufactured and evaluated to the highest standards.

#### **QUALITY ASSURANCE**

We do not just talk about quality, we are quality! Our pipette tips are ISO-certified and undergo rigorous quality standard testing. All our tips are independently certified free from RNase, DNase, pyrogens, Human DNA, ATP, and Endotoxin.

#### **CLEANROOM MANUFACTURING GUARANTEE**

Our LABBEE tips are manufactured in a 100,000 Class Clean Room facility. This ensures that they have passed the most stringent specifications for compatibility, cleanliness and consistency to meet the highest quality standards everywhere.

#### **WE DELIVER PRECISE RESULTS WITH ECONOMICAL PRICING, EVERY TIME**

Precision can only be delivered with the right products. Pipette tips are amongst the most commonly used products in a laboratory, so for anyone to settle for anything else other than top quality is never a choice.

The highest quality and economical pricing are not mutually exclusive and our LABBEE range ensures that you never have to make that compromise ever again.

#### **PRODUCT WHEN YOU NEED IT**

Being out of stock is unacceptable and leads to unnecessary disruptions. LABBEE ensures a constant stream of product to avoid such situations and to allow you to focus on what you do best, safe in the knowledge that we will always have the products you need available.



*Designed to meet the highest standards of precision and clarity*



## MANUAL TIPS

Item	SKU	Volume	Sterile	Conductive	Filter	Tip Color	Packaging	Number of units/case
Sterile Tips (manual tip)	LBFT10010	10 µL (31mm)	Yes	Yes	Yes	Clear	96/rack, 50 racks/case	4800
	LBFT10011	10 µL (39mm)	Yes	Yes	Yes	Clear	96/rack, 50 racks/case	4800
	LBFT10012	10 µL (46mm)	Yes	Yes	Yes	Clear	96/rack, 50 racks/case	4800
	LBFT10020	20 µL	Yes	Yes	Yes	Clear	96/rack, 50 racks/case	4800
	LBFT10100	100 µL	Yes	Yes	Yes	Clear	96/rack, 50 racks/case	4800
	LBFT10200	200 µL	Yes	Yes	Yes	Clear	96/rack, 50 racks/case	4800
	LBFT11000	1000 µL	Yes	Yes	Yes	Clear	96/rack, 50 racks/case	4800
	LBFT11250	1250 µL extra long	Yes	Yes	Yes	Clear	96/rack, 50 racks/case	4800

## ROBOTIC TIPS

Tecan Compatible Conductive Tips	LBTCBF0050	50 µL, conductive filter tips, sterile	Yes	Yes	Yes	Black	96/rack, 24 racks/case	2304
	LBTCBF0200	200 µL, conductive filter tips, sterile	Yes	Yes	Yes	Black	96/rack, 24 racks/case	2304
	LBTCBF1000	1000 µL, conductive filter tips, sterile	Yes	Yes	Yes	Black	96/rack, 24 racks/case	2304
Hamilton Compatible Conductive Tips	LBHTBF0050	50 µL, conductive filter tips, sterile	Yes	Yes	Yes	Black	96/rack, 24 racks/case	2304
	LBHTBF0300	300 µL, conductive filter tips, sterile	Yes	Yes	Yes	Black	96/rack, 24 racks/case	2304
	LBHTBF1000	1000 µL, conductive filter tips, sterile	Yes	Yes	Yes	Black	96/rack, 24 racks/case	2304





Victory Scientific  
PO Box 307 | Sewell, NJ 08080  
(856) 530-9990  
cs@victoryscientific.com

---

[victoryscientific.com](http://victoryscientific.com)