



50 x Disposable Face Mask

**All raw materials are made in USA
All manufacturing is made in USA**



PRODUCT DETAILS:

- **ASTM Level 1**
- **3-Ply**
- **Non-sterile**
- **Non-woven materials**
- **Outer layers are spun bond polypropylene**
- **Inner layer is melt-blown polypropylene**
- **Ear loop design**
- **Nose-piece wire forms to face to provide better seal**

Ordering Information

CATALOG #:

WFLMSK-1000

DESCRIPTION:

**ASTM Level 1, 3-Ply, Disposable, Protective Face Mask
with Earloops**

BREAKDOWN OF QUANTITIES:

50 Masks Per Box

32 Boxes Per Case

1600 Masks Per Case

PO Box 307

Sewell, NJ 08080

cs@victoryscientific.com

victoryscientific.com

(856) 530-9990

