



Mini Centrifuges

In Pink, Ivory-White and Blue

Our mini centrifuges are essential laboratory equipment used in clinical, environmental and industrial research. These personal centrifuges are ideal for quick spins, microfiltration of samples, cell separations and many other routine laboratory procedures.



Accessories



Standard Rotor for 6x1.5/2.0ml Tubes
Adaptors for 0.2ml Tubes (6 Pcs/Pack)

Rotor for 16x0.2ml Tubes
Adaptors for 0.5ml Tubes (6 Pcs/Pack)



Features of the Mini Centrifuges

- Ergonomically designed with a small footprint.
- High grade plastic cover for clear visibility.
- Smooth and quiet operation.
- Automatic stop when lid is opened.
- Digitally controlled, precise rotating speed.
- Fits a variety of tubes with included accessories.



Mini Centrifuges

Specifications

Name	Pink Mini Centrifuge	IvoryW Mini Centrifuge	Blue Mini Centrifuge
Model	MLX-104	MLX-106	MLX-110
Speed	4000 rpm	6000 rpm	10000 rpm
G-force (RCF)	900 xg	2000 xg	5000 xg
Noise Level	≤60 dB	≤60 dB	≤60 dB
Body Color	Pink	Ivory-White	Blue
Capacity	All : 6 x 1.5 ml / 6 x 0.5 ml / 16 x 0.2 ml		
Accessories Included	All : PCR Strip Rotor 16 x 0.2 ml & Standard Rotor 6 x 1.5 ml		
Ambient Temperature	5°C~40°C	5°C~40°C	5°C~40°C
Relative Humidity	≤80%	≤80%	≤80%
Power Supply	All : AC 100~240 V, 50 / 60 Hz		
Max Power	8W	10W	30W
Dimensions (WxDxH)	All : 155 x 140 x 115 mm (6.1 x 5.5 x 4.53")		
Net Weight	≤1 kg (2.2 lb)	≤1 kg (2.2 lb)	≤1 kg (2.2 lb)
Safety Feature	All : Stops When Lid is Opened		
Protection Level	IP21	IP21	IP21
Certifications	CE	CE	CE

Ordering Information

Cat. No.	Description	Pack
MLX-104	Pink Mini Centrifuge with Standard Accessories	1 SET
MLX-106	Ivory-White Mini Centrifuge with Standard Accessories	1 SET
MLX-110	Blue Mini Centrifuge with Standard Accessories	1 SET

Victory Scientific
PO Box 307
Sewell, New Jersey
(856) 530-9990
cs@victoryscientific.com
victoryscientific.com