



Low Profile Magnetic Stirrer

The magnetic stirrer is an essential instrument for every laboratory for mixing solution. It has a wide range of speed adjustment to meet the demands of different liquid mixing. It is designed with small size and high grade of waterproofing.



Accessories



Stirrer Bar Retriever: Ø7 x 300 mm

Stirrer Bar: Ø8 x 50 mm

Stirrer Bar: Ø10 x 40 mm



Stirrer Bar: Ø12 x 30 mm

Stirrer Bar: Ø7 x 20 mm

Stirrer Bar: Ø6 x 15 mm

Stirrer Bar: Ø12 x 50 mm



Features of the Low Profile Magnetic Stirrer

- Elegant design with high level of waterproof protection.
- Motorless driven mechanism for quiet and maintenance free operation.
- Sleek low-profile design.



Low Profile Magnetic Stirrer

Name	Low Profile Magnetic Stirrer
Model	MS-01FUA
Speed	15 ~ 1500 rpm
Auto. Reverse Rotation	Reversed Every 30 Seconds
Max. Length of Stir Bar	≤ 30 mm (1.2")
Platform Dimensions	110 x 110 mm (4.3 x 4.3")
Capacity	1 x 500 ml Bunsen Beaker or Flask
Max. Load	2 kg (4.4 lb)
Accessories Included	Stirrer Bar $\varnothing 7$ x 20 mm (0.28 x 0.79")
Ambient Temperature	5°C ~ 40°C
Relative Humidity	$\leq 80\%$
Power Supply	AC 100 ~ 240V, 50 / 60 Hz
Max. Power	6W
Platform Material	ABS
Dimensions (LxWxH)	175 x 125 x 20.5 mm (6.9 x 4.9 x 0.8")
Net Weight	0.3 kg (0.7 lb)
Protection Level	IP45
Certifications	CE, RoHS, WEEE

Ordering Information

Cat. No.	Description	Pack
MS-01FUA	Low Profile Magnetic Stirrer	1 SET

Victory Scientific
PO Box 307
Sewell, New Jersey
(856) 530-9990
cs@victoryscientific.com
victoryscientific.com